

PRESTO® A80t

Heating a 20 liters reactor from -60 °C to +20 °C

Objective

This case study tests the heating power of PRESTO® A80t with a 20 liters glass reactor. The PRESTO® A80t is connected to the reactor via two 1 m metal tubings. The PRESTO® A80t is programmed to heat up from -60 °C to +20 °C.



Environment

Room temperature +20 °C
 Humidity 45%
 Voltage 208 V / 60 Hz

Test Conditions


JULABO unit	PRESTO® A80t
Cooling power	+20 °C 1.2 kW 0 °C 1.2 kW -20 °C 1.1 kW
Heating capacity	3.4 kW
Band limit	with
Flow pressure	0.5 bar
Bath fluid	Thermal HL 80
Reactor	20 liters glass reactor (Chemglass) filled with 19 l Ethanol
Jacket volume	8 l
Control	External (ICC)



Test Results

The PRESTO® A80t heating process from -60 °C to +20°C in 1 h 30 min without overshoot.

Measured with *EasyTEMP Professional*



- Setpoint
- Temperature in reactor's interior
- Temperature in reactor's jacket

Tip

Take advantage of our wide range of accessories. The M+R adapter enables you to display and record an additional temperature.



Tip

Use the free of charge *EasyTEMP* software to control the units with the PC and to show the temperature curves graphically.

EasyTEMP

